

## COVID-19 Home Test Kit Important Notice

Thank you for taking the home test for coronavirus and helping to protect your community.

You will find the full instruction inside the test kit

Please ensure:

1. You register your test as per the instructions below
2. You place the bar code vertically on the plastic vial that corresponds to the barcode on the Biohazard Bag
3. The plastic vial is placed in the clear zip-lock bag **with** the absorbent pad
4. The clear zip-lock bag is placed inside the Biohazard Bag
5. The Biohazard Bag is sealed using the silver seal
6. Assemble the cardboard box provided and place the Biohazard Bag securely inside
7. Apply the Security Seal to secure the bag closed

### IMPORTANT

To get your test you **MUST**:

Register outline at

[www.test-for-coronavirus.service.gov.uk/register-kit](http://www.test-for-coronavirus.service.gov.uk/register-kit) or call 199

**Our staff will return today to collect your completed tests  
If for any reason your test is not collected today you can post  
your test kit in your nearest priority postbox or have it collected  
by a courier.**

Please see page 17 of the 'COVID-19 Home Test Step-by-step guide' for details postal instructions